


Harnessing Photoplethysmography and Artificial Intelligence for Early Diagnosis of Congenital Heart Disease in Newborns

Exploring an innovative approach to combine non-invasive photoplethysmography (PPG) technology with advanced artificial intelligence algorithms to improve the early detection of congenital heart disease (CHD) in newborns



In The Name of God

Harnessing Photoplethysmography and Artificial Intelligence for Early Diagnosis of Congenital Heart Disease in Newborns

Seyyed Ali Najafi

Summary

The integration of photoplethysmography (PPG) and artificial intelligence (AI) presents a transformative potential for diagnosing congenital heart disease (CHD) in newborns. PPG is a non-invasive optical technique that measures blood volume changes in the microvascular bed of the skin, providing real-time data on cardiovascular health. This technology is particularly advantageous for neonatal care, where traditional diagnostic methods can be invasive or logistically challenging. By combining PPG with AI algorithms, healthcare professionals can enhance the accuracy and speed of CHD diagnosis, facilitating earlier interventions that are critical for improving patient outcomes.

The notable intersection of PPG and AI is driven by advances in machine learning and deep learning techniques, which allow for the sophisticated analysis of complex physiological data. Recent studies have shown that AI models can achieve high accuracy rates in identifying cardiovascular issues from PPG signals, with some models reporting detection exceeding 96% accuracy rate.

This capability is particularly valuable in resource-limited settings, where access to advanced imaging technologies may be restricted, and mobile health applications can enable continuous monitoring of newborns at risk for CHD.

As research continues to evolve, the collaboration among clinicians, AI researchers, and bioinformaticians will be essential in overcoming existing barriers and maximizing the benefits of PPG and AI for diagnosing CHD in newborns. Future advancements in signal processing, machine learning methodologies, and multidisciplinary partnerships hold the promise of significantly improving early detection and management of congenital heart defects, ultimately enhancing the quality of care for vulnerable pediatric populations.



Congenital Heart Disease (CHD)

Importance of Early Detection

The birth prevalence of congenital heart disease globally increased to a maximum of 9.410/1000 in 2010-17, with Africa reporting the lowest prevalence and Asia the highest. Early detection of congenital heart disease (CHD) is crucial for improving outcomes in affected infants. Early detection of congenital heart disease in newborns is important, as treatment outcomes are related to the time of diagnosis. Early diagnosis of congenital heart diseases contributes to reducing the mortality of newborns and improving therapeutic plans.

The implementation of screening protocols, particularly through methods such as photo plethysmography, has shown promise in identifying critical congenital heart defects (CCHDs) early in life, potentially allowing for timely intervention and management.

Common Types of Congenital Heart Defects

There are 18 recognized types of congenital heart defects, each affecting the heart's anatomy and function in different ways.

- Atrial Septal Defect (ASD): This condition involves an abnormal opening between the two upper chambers of the heart, which can lead to excessive blood flow to the lungs.
- Ventricular Septal Defect (VSD): Characterized by a hole in the ventricular septum, this defect allows blood to flow from the left ventricle back into the right ventricle, increasing pulmonary blood flow and pressure.
- Coarctation of the Aorta (CoA): This defect involves the narrowing of the aorta, which can obstruct blood flow to the lower body and raise blood pressure above the constriction.
- Hypoplastic Left Heart Syndrome (HLHS): A complex condition involving several abnormalities that affect normal blood flow through the heart, leading to inadequate blood supply to the body.
- Tetralogy of Fallot (TOF): This condition includes four defects that result in insufficient blood flow to the lungs and low oxygen levels in the bloodstream.
- Transposition of the Great Arteries: A serious defect in which the positions of the aorta and pulmonary artery are switched, leading to inadequate oxygenation of blood.



Image1. Critical conditions like Tetralogy of Fallot (TOF), Hypoplastic Left Heart Syndrome (HLHS), and Transposition of the Great Arteries highlight the urgent need for early detection and intervention to improve survival and outcomes in affected newborns

Diagnostic Accuracy

The accuracy of diagnoses in congenital heart disease (CHD) is critical for effective management, yet it poses significant challenges. One key method used in diagnosis, transthoracic echocardiography (TTE), has been noted to be difficult to perform and interpret due to factors such as the complexity of cardiac anatomy and the varying degrees of cooperation from children.

Transthoracic echocardiography (TTE) is an invaluable diagnostic tool for congenital heart disease (CHD), but it has notable limitations. TTE is not routinely performed for all children and is reserved for those with clinical signs or symptoms of CHD. This selective use may delay the detection of subtle or asymptomatic cases.

Even when performed, TTE may misinterpret complex cardiac anatomies, leading to delayed diagnoses or unnecessary invasive procedures. Additionally, it fails to provide sufficient diagnostic information in approximately 27% of cases, requiring supplementary imaging such as cardiac catheterization or CT angiography. These limitations highlight the need for complementary tools and techniques to enhance diagnostic accuracy in CHD.



Image2. Finger-based SpO₂ devices showcase the practical, day-to-day application of PPG in monitoring vital signs.

Photoplethysmography (PPG) Technology

Photoplethysmography (PPG) is a non-invasive optical technique used to measure blood volume changes in the microvascular bed of the skin. This method is based on the principles of light absorption, scattering, and transmission properties of biological tissues when illuminated by a specific wavelength of light. The term PPG combines "photo," meaning light; "plethysmo," referring to volume; and "graphy," meaning recording. The foundational work of Hertzman in 1937 highlighted the relationship between cardiac activity and changes in light detected via backscattering, establishing PPG as a viable technique for monitoring blood volume fluctuations in targeted areas.

Measurement Principles

The PPG technique involves a light source and a photodetector placed on the skin. Light emitted from the source penetrates the skin and is absorbed by the blood, with the intensity of the reflected light captured by the detector. This process results in a waveform that reflects the pulsatile nature of blood flow corresponding to the cardiac cycle. During the diastolic phase, blood volume, arterial diameter, and hemoglobin concentration are at their minimum, maximizing light detection by the photodetector, whereas during the systolic phase, the opposite occurs, resulting in minimal light detection. This creates a time-series waveform that represents the changes in blood volume corresponding to each heartbeat. The resulting PPG data includes information on the amplitude and frequency of these fluctuations, which reflects both the heart rate and the variability in blood circulation. The data is typically recorded at a high sampling rate to capture detailed cardiac cycles and can be used to monitor physiological parameters such as heart rate, oxygen saturation, and blood flow dynamics.




Image3. PPG technology generates data on oxygen saturation (SpO_2) by analyzing blood volume changes in response to light absorption

Applications in Healthcare

PPG technology has gained popularity in both clinical and mobile health applications due to its non-invasive nature and ease of use. It is increasingly utilized for monitoring various physiological parameters, such as heart rate, respiratory rate, and stress levels. PPG is increasingly utilized in neonatal care for non-invasive monitoring of vital signs, particularly heart rate (HR) and respiratory rate (RR). This technology offers significant advantages in clinical settings, especially in neonatal intensive care units (NICUs), where traditional monitoring methods can be intrusive. The versatility of PPG allows for its implementation in wearable devices that continuously assess users' physiological states, thus enhancing daily health monitoring and clinical diagnostics. Additionally, the integration of PPG with artificial intelligence (AI) is emerging as a promising area for diagnosing conditions such as congenital heart disease (CHD) in newborns, potentially improving early detection and intervention.

Signal Processing Challenges

Despite its advantages, PPG signals are susceptible to noise and artifacts, particularly from movement and environmental factors. The quality of PPG signals can degrade due to low perfusion, ambient light interference, and individual differences such as skin type. Addressing these challenges is crucial for enhancing the reliability of PPG in clinical settings. Current research emphasizes the need for improved preprocessing techniques, including machine learning-based approaches, to differentiate between analyzable and non-analyzable segments of PPG data. As the field evolves, there is a continuous demand for refining PPG technology to ensure its effectiveness in both clinical and mobile health environments.

Artificial Intelligence in Healthcare

Artificial intelligence (AI), particularly through the use of machine learning (ML) and deep learning (DL) techniques, has shown significant potential to transform healthcare, particularly in the diagnosis and treatment of congenital heart disease (CHD) in newborns. AI systems are capable of analyzing vast amounts of complex medical data, including imaging, genetic information, and clinical records, enabling them to identify patterns and make predictions that may exceed human capabilities.




Image4. Ai is Transforming medical care through advanced algorithms that enhance diagnosis, treatment, and patient outcomes

Predictive Analytics and Early Intervention

One of the key applications of AI in healthcare is predictive analytics. ML models can analyze patient records and risk factors to predict the likelihood of CHD, allowing for earlier interventions. This capability is especially crucial in overburdened healthcare systems, where prioritizing high-risk patients for treatment can improve outcomes.

Enhancing Diagnostic Accuracy

AI has advanced the accuracy of CHD diagnosis through sophisticated algorithms that analyze imaging data. For instance, deep learning methods have been utilized to enhance imaging techniques and assist in surgical planning. Research indicates that various AI models, including support vector machines (SVM) and neural networks, have significantly improved detection performance in terms of accuracy, sensitivity, and specificity when analyzing multimodal features related to CHD.

Combining PPG and AI for CHD Diagnosis

The integration of Photoplethysmography (PPG) and artificial intelligence (AI) presents a promising avenue for improving the diagnosis of congenital heart disease (CHD) in newborns. PPG is a non-invasive optical technique that can be used to detect blood volume changes in the microvascular bed of tissue, making it a valuable tool for monitoring cardiovascular conditions.

Advantages of PPG in CHD Diagnosis

PPG's non-invasive nature allows for continuous monitoring of infants, which is particularly advantageous in a neonatal setting where traditional methods, such as echocardiograms, may be challenging to perform repeatedly due to their complexity and need for skilled personnel. The ability to collect data in real-time can lead to earlier detection of potential heart defects, allowing for timely intervention and management.

AI-Driven Analysis of PPG Data

Recent studies have demonstrated that artificial intelligence (AI) can significantly enhance the interpretation of PPG signals for diagnosing congenital heart disease (CHD). AI techniques are capable of identifying subtle patterns in PPG data that traditional methods may overlook. For instance, models trained on large datasets of newborns have achieved remarkable accuracy rates, with some exceeding 96% in detecting cardiovascular abnormalities. These results are backed by rigorous validation processes using diverse patient populations. AI-powered analysis also allows for rapid, consistent, and non-invasive screening, which is particularly valuable in neonatal care. By integrating AI with PPG technology, we can offer a scalable solution that improves diagnostic precision and ensures timely intervention for at-risk newborns, even in resource-limited settings.





Image5. SpO₂ measurements from the foot of a newborn provide crucial data for developing and training AI algorithms

Data Collection and Required Data Types

For the successful integration of AI with PPG signals to diagnose congenital heart disease (CHD), a diverse dataset is essential. The types of data we will gather include:

1. PPG Signals:

- Continuous, real-time PPG data collected from newborns using a wearable sensor.
- PPG data must include time-series recordings of blood volume changes in response to cardiac cycles, typically captured at a sampling rate of 100 Hz.
- Signal variations from healthy infants and those diagnosed with CHD will be recorded to help the model learn the distinct patterns associated with heart abnormalities.

2. Clinical and Demographic Information:

- Health status: Details on the newborn's condition, including any diagnosed heart defects or abnormalities (e.g., Ventricular Septal Defect, Tetralogy of Fallot).
- Medical History: Family history of congenital heart disease, gestational age, birth weight, and any previous medical interventions.
- Physical Exam Data: Recorded physical findings such as heart rate, blood pressure, and oxygen saturation levels.

3. Additional Physiological Signals:

- In some cases, complementary data from ECG or pulse oximetry will be integrated with PPG signals to enhance the diagnostic model's performance. These signals will help in cross-validating the results derived from PPG analysis.

4. Labeling and Diagnosis:

- Each dataset will be labeled with diagnostic outcomes (e.g., presence or absence of CHD, type of defect). This will involve collaboration with neonatologists and cardiologists to ensure accurate, clinically-relevant labeling.
- The final dataset will be divided into two categories: infants diagnosed with specific CHD types and healthy controls for training and validation.




Image6. Finger-mounted SpO2 devices enable non-invasive and continuous monitoring in newborns

Data Preparation and Model Training

The training process for these AI models involves several critical steps, including data standardization, feature extraction, and the implementation of advanced training methodologies. Data collected from PPG devices must first undergo cleaning and filtering to eliminate noise and irrelevant variables. This ensures that the models trained on this data can accurately learn the patterns associated with CHD. Additionally, leveraging large datasets of labeled PPG recordings from both healthy infants and those diagnosed with CHD enhances the model's ability to generalize and perform accurately in clinical settings.

Challenges

Despite the advantages, challenges remain in the widespread adoption of PPG and AI for CHD diagnosis. Issues such as variations in PPG signal quality due to external factors (e.g., motion artifacts) and the need for robust validation of AI models in diverse populations must be addressed. Future research should focus on refining these AI algorithms, enhancing their robustness, and integrating them into clinical workflows to facilitate their use by healthcare professionals. Additionally, exploring the combination of PPG data with other physiological signals may further improve diagnostic accuracy and patient outcomes in CHD detection.

